

MIB Trousse de chirurgie pour Mini Implant complète

Pilier de cicatrisation



10 N.cm	TOURNEVIS MANUEL	MANDRIN C.A
	▶ TMH12-C	▶ MTH12-C
	▶ TMH12-L	▶ MTH12-L

Prothèse Vissée



Prothèse Scellée



15 N.cm	TOURNEVIS MANUEL	MANDRIN C.A
	▶ TMH12-C	▶ MTH12-C
	▶ TMH12-L	▶ MTH12-L

15 N.cm	TOURNEVIS DYNAMIC C.A
	▶ L:18MM 43.618.201.01-2
	▶ L:24MM 43.624.201.01-2
	▶ L:32MM 43.632.201.01-2

Amovible



10 N.cm	TOURNEVIS MANUEL	MANDRIN C.A
	▶ TMH12-C	▶ MTH12-C
	▶ TMH12-L	▶ MTH12-L


Transferts



Scanbody



Analogues



Vis de laboratoire

